


**Quick and Easy Calibration** - of Installed DAP (KAP) and Patient Dose Systems

*What's it all about?*



Foam Elevation Support Stand



*Patient Dose Control*



**TRACEABLE MEASUREMENTS** - The PDC is a reference class instrument for "field calibration" of patient dose measurement and control systems thus ensuring the validity of inter-institution patient dose comparisons.

**FLEXIBLE AND CONVENIENT** - Use the PDC with a phantom to simulate patient imaging conditions including scattered radiation or mount the PDC on its support stand to measure air Kerma.

**SIMPLE TO USE** - The PDC displays DAP (KAP) and dose rate during an exposure then automatically switches to display accumulated DAP (KAP) and dose on exposure completion.

**DEPENDABLE** - A tough ABS plastic housing protects the ion chambers and electronics that incorporate several patented features to ensure long term stability.

**DAP CALIBRATION TEMPLATE**

mA	Rm. 1 Dose (mGy)	Rm. 1 DAP (µGy·m²)	Sample	PDC Dose (mGy)	PDC DAP (µGy·m²)	% DOSE VARIATION	% DAP VARIATION
50	0.3	6	1	0.30	6.58	1%	-9%
100	0.6	12	2	0.59	12.69	2%	-5%
200	1.1	25	3	1.17	25.88	-6%	-3%
400	2.3	50	4	2.30	51.81	0%	-3%

**THE GOLD STANDARD IN RADIATION MEASUREMENT**

# PDC KEY FEATURES AND BENEFITS:

## KEY FEATURES

**Complete DAP meter assessment:**

**Symmetrical Response:**

**Patient Dose System assessment:**

1. Dose/DAP Air Kerma Calibration:
2. Entrance (skin) Dose / DAP Calibration and QA:
3. Image Receptor Dose/ DAP Calibration and QA:

**Optional Remote Control Software:**

**Optical and radiographic alignment markers:**

## BENEFITS

Measures DAP, DAP-rate over a full range of field sizes and beam qualities

Can be used with under-couch tubes without the need for inversion

Measure Dose and Dose-rate

Measure DAP and dose linearity with dose and field size


Surface dose calibration with phantom and establish reference DAP and dose levels

Verify Image Receptor dose measurements and AEC linearity testing

Automatic data capture with customizable templates


Setting reference field sizes made simple

## APPLICATIONS




### DOSE / DAP AIR KERMA CALIBRATION

- DAP and dose linearity with dose and field size.
- Sensitivity testing.
- DAP and dose calibration at the patient plane or at a reference distance.



### ENTRANCE (SKIN) DOSE / DAP CALIBRATION AND QA

- Entrance surface dose calibration for different examinations using a phantom.
- Establishing examination related reference DAP and dose levels.



### IMAGE RECEPTOR DOSE / DAP CALIBRATION AND QA

- Image receptor dose measurement.
- AEC linearity testing.

## SPECIFICATIONS / TECHNICAL DATA:

### Display Range

Dose area product  
Air kerma

(0.01... 99 999 999)  $\mu\text{Gy}\cdot\text{m}^2$   
(0.001... 99 999 999) mGy

### Accuracy

DAP and Air kerma

Inclusive of all uncertainties (temperature, pressure, rate, area and beam quality)  $\pm 10\%$   
Under reference conditions (10 mGy/min, 15 X 15 cm field, 80 kVp, 2.5 mm Al filtration)  $\pm 7.5\%$

### Digital Resolution

Dose area product  
Dose area product rate  
Air kerma  
Air kerma rate

0.01  $\mu\text{Gy}\cdot\text{m}^2$   
1  $\mu\text{Gy}\cdot\text{m}^2/\text{min}$   
0.001 mGy  
0.1 mGy/min

### Rated range of use

Tube voltage  
Dose area product:  
(low rate range)  
(high rate range)  
Air kerma rate

(40 ... 150) kV  
  
(1 ...  $1\cdot 10^4$ )  $\mu\text{Gy}\cdot\text{m}^2/\text{min}$   
( $2\cdot 10^3$  ...  $9\cdot 10^5$ )  $\mu\text{Gy}\cdot\text{m}^2/\text{min}$   
(0.2 ...  $9\cdot 10^3$ ) mGy/min (at the position of the chamber)

### Automatic Temperature and Pressure Correction

Pressure  
Temperature  
Air humidity

(80.0 ... 106.0) kPa  
(+10 ... +40) °C  
(10 ... 80) % rel. humidity (max. 20 g/m<sup>3</sup>)

### Ionization chamber

Response versus radiation quality (50kV ... 150 kV, norm. to 100kV; acc. IEC 60580)  
Quality equivalent filtration (70kV)  
Active area  
Dose area product  
Air kerma

$\pm 3\%$   
0.6 mm Al  
max. (300 x 300) mm<sup>2</sup>  
min. (100 x 100) mm<sup>2</sup>

### Power supply

Internal rechargeable battery pack  
operation time (state of charge: 100%)

Li-ION, 2 cells  
> 8 h

### Serial Interface

Protection class (acc. IEC 60529)

Weight  
Dimension

USB  
IP 41  
2.32 kg  
350 mm x 410 mm x 35 mm  
(L x W x H)

All specifications subject to change.



### PORTABLE CONVENIENT CARRY CASE

Foam Elevation Support Stand fits inside the interior of the Carry Case lid.

