

PRONE HEAD AND NECK IMMOBILIZATION SYSTEM



- ▶ Adjustable from positive 20° to negative 20° in 5° increments
- ▶ Supine adapter elevates head and provides tilt ideal for patients that cannot lay flat
- ▶ Compatible with all CDR Freedom Immobilization



Supine adapter shown ID# PH08-SH



PH08
Prone Head and Neck Immobilization System

PRONE FACE MASK MAKER



- ▶ Used as the initial setup for making the Prone face mask when setting up patient on PH08 system
- ▶ Quick adjust push button adjustment for setting face mask position



PH08-SD
Hands-free Prone Mask Maker for Prone Head and Neck Immobilization System