

BEEKLEY MR SPOTS

MR-SPOTS:

Beekley MR-SPOTS® **image on all MRI sequences** so that areas of interest are clearly identified on the scan.

MR-SPOTS® tubes provide a sharp, distinct image on film while MR-SPOTS packets help reduce marker indentation.




3.0cm Radiance® filled tube



MR-SPOTS®
order code: 122
contents per box: 40
1.5cm Radiance® filled tube
perfect for extremities, soft tissue masses, and breast MRI



MR-SPOTS®
order code: 121
contents per box: 40
3.0cm Radiance filled tube
perfect for thoracic spines and larger areas of interest



1.75cm Radiance filled packet



MR-SPOTS®
order code: 185
contents per box: 40
1.75cm Radiance filled packet
perfect for extremities, soft tissue masses, and breast MRI
designed for optimum patient comfort



MR-SPOTS®
order code: 184
contents per box: 40
5.0cm Radiance filled packet
perfect for thoracic spines and larger areas of interest
designed for optimum patient comfort