

## WIPE TEST KITS



**W**ipe Test Kits may be used to perform wipe tests as indicated in NRC and Agreement State Regulations. They provide an efficient, convenient means of sampling contaminated areas with radioactivity on either

wet or dry surfaces.

Each box of wipes includes 500 record folders which may be used to: (a) identify each sample wipe, (b) prevent cross-contamination of smears, (c) transfer the wipe to the counting facility, and (d) store the wipe until discarded.

Wipes are 1.75" dia cloth discs packaged in a box of five hundred.

**006-350** Wipe Test Kits, 1.75" dia, 500/pkg